


Product Datasheet

Human Rantes ELISA Kit (orb50098)

Description	Human Rantes ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S24-S91
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
MW	9990 MW
Uniprot ID	D0EI67
Sensitivity	3 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum is 6.8-17.8 ng/ml.
Expiration Date	Please enquire.



Human
Rantes
PicoKine
ELISA Kit
standard...