



Product Datasheet

Neutrophil Elastase antibody (orb500993)


Description	Neutrophil Elastase antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Neutrophil Elastase/ELANE (101-200/267 aa)
Target	Neutrophil Elastase
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26 kDa
Uniprot ID	P08246
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Western blot analysis of Spinal cord (Mo...)



Western blot analysis of Lane 1: Bone (M...)



Western blot analysis of HL60(Human) Cel...