



Product Datasheet

Aquaporin 4 antibody (orb501006)


Description	Aquaporin 4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human AQP4 (271-323/323 aa)
Target	Aquaporin 4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	P55087
Dilution Range	IHC-P=1:100-500, Flow-Cyt=1µg /test, WB=1:500-2000, IHC-F=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of Muscle(Mouse) l...



Western blot analysis of mouse lung lysa...



Immunohistochemical staining of human la...