


Product Datasheet

Human P-Selectin / CD62P ELISA Kit (orb50101)

Description	Human P-Selectin / CD62P ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: W42-A771
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	90834 MW
Uniprot ID	P16109
Sensitivity	5 pg/ml
Dilution Range	Dilution ratio of 1:20, concentration in serum and plasma is 51-113 ng/ml.
Expiration Date	Please enquire.



Human P-
Selectin
PicoKine
ELISA Kit
stan...