


Product Datasheet

Rat TGF Beta 1 ELISA Kit (orb50104)

Description	Rat TGF Beta 1 ELISA Kit
Reactivity	Rat
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: A279-S390
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (EDTA) and urine.
MW	44329 MW
Uniprot ID	P17246
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:500, concentration in serum and plasma is 24-85 ng/ml.
Expiration Date	Please enquire.



Rat TGF beta
1 PicoKine
ELISA Kit
standa...