


Product Datasheet

Human TIE2/TEK ELISA Kit (orb50107)

Description	Human TIE2/TEK ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A23-K745
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	125830 MW
Uniprot ID	Q02763
Sensitivity	5 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is 14-90 ng/ml.
Expiration Date	Please enquire.



Human TIE2
PicoKine
ELISA Kit
standard c...