


Product Datasheet

Human TIMP-1 ELISA Kit (orb50108)

Description	Human TIMP-1 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C24-A207
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and saliva.
MW	23171 MW
Uniprot ID	P01033
Sensitivity	5 pg/ml
Dilution Range	Dilution ratio of 1:1000, concentration in serum is 71.8-137 ng/ml.
Expiration Date	Please enquire.



Human TIMP-1 PicoKine ELISA Kit standard...