


Product Datasheet

Human TIMP-2 ELISA Kit (orb50109)

Description	Human TIMP-2 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: C27-P220
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and saliva.
MW	24399 MW
Uniprot ID	P16035
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is 18-393 ng/ml.
Expiration Date	Please enquire.



Human TIMP-2 PicoKine ELISA Kit standard...