

Product Datasheet

PRKD1 (phospho-Ser738) antibody (Biotin) [Out of stock] (orb501125)

Description

Rabbit polyclonal antibody against

conjugated to Biotin

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Biotin

Tested

IHC-P

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human PRKD1 around the phosphorylation site of Ser738 EK(p-S)FR

Preservatives

0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

102kDa

Uniprot ID
Q15139
Dilution Range

WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.