

Product Datasheet

Mouse TNF Alpha/Tumor necrosis factor ELISA Kit (orb50113)

Description Mouse TNF Alpha/Tumor necrosis factor ELISA Kit

Reactivity Mouse

Range 15.6 pg/ml - 1,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: E.coli; Immunogen sequence: L80-L235

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only

Sample Types cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).


MW 25896 MW

Uniprot ID **P06804**

Sensitivity 1 pg/ml

Dilution Range Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 15.6 pg/ml.

Expiration Date Please enquire.



Mouse
TNF alpha
ELISA Kit