


Product Datasheet

Human sTNFsR II ELISA Kit (orb50114)

Description	Human sTNFsR II ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L23-D257
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA, citrate) and urine.
MW	48291 MW
Uniprot ID	P20333
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:20, concentration in serum and plasma is 802-3804 pg/ml.
Expiration Date	Please enquire.



Human
sTNFsR II
PicoKine
ELISA Kit
stand...