


Product Datasheet

Human uPA / PLAU / URK ELISA Kit (orb50115)

Description	Human uPA / PLAU / URK ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: M1-L431
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	48507 MW
Uniprot ID	P00749
Sensitivity	5 pg/ml
Dilution Range	For serum, dilution ratio is 1:1 and concentration is 500-1700 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 1000-2000 pg/ml.
Expiration Date	Please enquire.



Human
uPA/PLAU
PicoKine
ELISA Kit
standa...