


Product Datasheet

Human VCAM-1 ELISA Kit (orb50117)

Description	Human VCAM-1 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: F25-E698
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	81276 MW
Uniprot ID	P19320
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:300, concentration in serum and plasma is around 500 ng/ml.
Expiration Date	Please enquire.



Human
VCAM-1
PicoKine
ELISA Kit
standard...