


Product Datasheet

Human VEGFR1/sFLT1 ELISA Kit (orb50122)

Description	Human VEGFR1/sFLT1 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S27-H687
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (EDTA).
MW	150769 MW
Uniprot ID	P17948
Sensitivity	30 pg/ml
Dilution Range	For serum, dilution ratio is 1:1 and concentration is 3.6-5 ng/ml. For plasma, dilution ratio is 1:1 and concentration is less than the lowest standard, 156 pg/ml.
Expiration Date	Please enquire.



Human VEGFR1/sFLT1 PicoKine ELISA Kit st...