


Product Datasheet

Human ACE/Cd143 ELISA Kit (orb50123)

Description	Human ACE/Cd143 ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L30-L1261
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin) and saliva.
MW	53750 MW
Uniprot ID	P12821
Sensitivity	5 pg/ml
Dilution Range	Dilution ratio of 1:20, concentration in serum and plasma is 29-253 ng/ml.
Expiration Date	Please enquire.



Human ACE
PicoKine
ELISA Kit
standard
cu...