

Product Datasheet

Human E-Cadherin/CDH1/Cadherin-1 ELISA Kit (orb50124)

Description Human E-Cadherin/CDH1/Cadherin-1 ELISA Kit

Reactivity Human

Range 156 pg/ml - 10,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: NS0; Immunogen sequence: D155-I707

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only

Sample Types cell culture supernatants, serum, plasma (heparin) , saliva and urine.


MW 97456 MW

Uniprot ID **P12830**

Sensitivity 15 pg/ml

Dilution Range Dilution ratio of 1:10, concentration in serum and plasma is 17-71 ng/ml.

Expiration Date Please enquire.



Human E-Cadherin PicoKine ELISA Kit stan...