

Product Datasheet

Cdc6 (phospho-Ser54) antibody (Biotin) (orb501254)

Description	Cdc6 (phospho-Ser54) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, ICC, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Cdc6 around the phosphorylation site of Ser54
Target	phospho-Cdc6 (Ser54)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q99741
Dilution Range	WB=1:500-2000,ICC=1:100,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.