


Product Datasheet

Rat TIMP-1 ELISA Kit (orb50129)

Description	Rat TIMP-1 ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C24-A217
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	23794 MW
Uniprot ID	P30120
Sensitivity	3 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is 5920-49680 pg/ml.
Expiration Date	Please enquire.



Rat TIMP-1
PicoKine
ELISA Kit
standard c...