


Product Datasheet

Human TNF Beta/Lymphotoxin-alpha ELISA Kit (orb50130)

Description	Human TNF Beta/Lymphotoxin-alpha ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: L35-L205
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	22297 MW
Uniprot ID	P01374
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Human TNF
beta
PicoKine
ELISA Kit
standa...