


Product Datasheet

Human VEGF-C ELISA Kit (orb50131)

Description	Human VEGF-C ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T103-R227
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	46883 MW
Uniprot ID	P49767
Sensitivity	3 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 500-1250 pg/ml.
Expiration Date	Please enquire.



Human
VEGF-C
PicoKine
ELISA Kit
standard...