


Product Datasheet

Human Adiponectin ELISA Kit (orb50133)

Description	Human Adiponectin ELISA Kit
Reactivity	Human
Range	1.56 ng/ml - 100 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E19-N244
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	26414 MW
Uniprot ID	A8K660
Sensitivity	60 pg/ml
Dilution Range	Dilution ratio of 1:5000, concentration in serum and plasma is around 90 ug/ml.
Expiration Date	Please enquire.



Human
Adiponectin
PicoKine
ELISA Kit
sta...