


Product Datasheet

Mouse Adiponectin ELISA Kit (orb50134)

Description	Mouse Adiponectin ELISA Kit
Reactivity	Mouse
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E18-N244
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	26809 MW
Uniprot ID	Q60994
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:2000, concentration in serum and plasma is 4.8-6.1 ug/ml.
Expiration Date	Please enquire.



Mouse
Adiponectin
PicoKine
ELISA Kit
sta...