


Product Datasheet

Human TSP2/THBS2/Thrombospondin-2 ELISA Kit (orb50136)

Description	Human TSP2/THBS2/Thrombospondin-2 ELISA Kit
Reactivity	Human
Range	625 pg/ml - 40,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: G19-I1172
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and human milk.
MW	129991 MW
Uniprot ID	P35442
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:10, concentration in serum and plasma is around 50 ng/ml.
Expiration Date	Please enquire.



Human
TSP2
PicoKine
ELISA Kit
standard c...