


Product Datasheet

Human Endoglin/CD105/ENG ELISA Kit (orb50137)

Description	Human Endoglin/CD105/ENG ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: E26-G586
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	70578 MW
Uniprot ID	P17813
Sensitivity	15 pg/ml
Dilution Range	Dilution ratio of 1:4, concentration in serum and plasma is 2-8 ng/ml.
Expiration Date	Please enquire.



Human
CD105
PicoKine
ELISA Kit
standard ...