


Product Datasheet

Human P-Cadherin-3 CDH3 ELISA Kit (orb50138)

Description	Human P-Cadherin-3 CDH3 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D108-G654
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	91418 MW
Uniprot ID	P22223
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum is around 8 ng/ml.
Expiration Date	Please enquire.



Human P-Cadherin PicoKine ELISA Kit stan...