


Product Datasheet

Human Cystatin C/CST3 ELISA Kit (orb50140)

Description	Human Cystatin C/CST3 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: M1-A146
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, urine and human milk.
MW	15799 MW
Uniprot ID	P01034
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:500, concentration in serum and plasma is around 1500 ng/ml.
Expiration Date	Please enquire.



Human Cystatin C
PicoKine
ELISA Kit
stan...