

## Product Datasheet

### NFKB1 (phospho-Ser927) antibody (Biotin) (orb501424)

|                            |                                                                                                            |
|----------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | NFKB1 (phospho-Ser927) antibody (Biotin)                                                                   |
| <b>Species/Host</b>        | Rabbit                                                                                                     |
| <b>Reactivity</b>          | Rat                                                                                                        |
| <b>Conjugation</b>         | Biotin                                                                                                     |
| <b>Tested Applications</b> | ELISA, IHC-Fr, IHC-P                                                                                       |
| <b>Immunogen</b>           | KLH conjugated Synthesised phosphopeptide derived from rat NFKB1 around the phosphorylation site of Ser927 |
| <b>Target</b>              | Phospho-NFKB1 (Ser927)                                                                                     |
| <b>Preservatives</b>       | Preservative: 15mM Sodium Azide,<br>Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4                           |
| <b>Form/Appearance</b>     | Liquid                                                                                                     |
| <b>Concentration</b>       | 1mg/ml                                                                                                     |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.                           |
| <b>Note</b>                | For research use only                                                                                      |
| <b>Isotype</b>             | IgG                                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                                 |
| <b>Uniprot ID</b>          | <b>P19838</b>                                                                                              |
| <b>Dilution Range</b>      | IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                            |