


Product Datasheet

Human CXCL9/Mig ELISA Kit (orb50143)

Description	Human CXCL9/Mig ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: T23-T125
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	14019 MW
Uniprot ID	Q07325
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 200 pg/ml.
Expiration Date	Please enquire.



Human
CXCL9
PicoKine
ELISA Kit
standard ...