


Product Datasheet

Human Lipocalin-2/NGAL ELISA Kit (orb50149)

Description	Human Lipocalin-2/NGAL ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q21-G198
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin) , saliva and urine.
MW	22588 MW
Uniprot ID	P80188
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is 40-100 ng/ml.
Expiration Date	Please enquire.



Standard Curve of Human NGAL ELISA Kit.