


Product Datasheet

Mouse Lipocalin-2/NGAL ELISA Kit (orb50150)

Description	Mouse Lipocalin-2/NGAL ELISA Kit
Reactivity	Mouse
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q21-N200
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	22875 MW
Uniprot ID	P11672
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:500, concentration in serum and plasma is 160-800 ng/ml.
Expiration Date	Please enquire.



Mouse
Lipocalin-
2/NGAL
PicoKine
ELISA Ki...