


Product Datasheet

Rat Lipocalin-2/NGAL ELISA Kit (orb50151)

Description	Rat Lipocalin-2/NGAL ELISA Kit
Reactivity	Rat
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q21-N198
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin) and urine.
MW	22476 MW
Uniprot ID	P30152
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Rat
Lipocalin-
2/NGAL
PicoKine
ELISA Kit ...