

Product Datasheet

Connexin 43 (phospho-Ser368) antibody (Biotin)
[Out of stock] (orb501518)

Description

Rabbit polyclonal antibody

Species/Host

Rabbit

Conjugation

Biotin

Tested

ELISA, IF, IHC-Fr, IHC-P

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human Connexin 43 around the phosphorylation site of Ser368

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

43kDa

Uniprot ID

P23242, Q6TYA7, Q6S9C0, P18246, P08050, P17302

Dilution Range

ELISA: 1:100-1000, IHC-P: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.