


Product Datasheet

Human IL-18/IL-1F4/Interleukin-18 ELISA Kit (orb50153)

Description	Human IL-18/IL-1F4/Interleukin-18 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Y33-D197
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	22326 MW
Uniprot ID	Q14116
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 67-429 pg/ml.
Expiration Date	Please enquire.



Human IL-18 PicoKine ELISA Kit standard ...