


Product Datasheet

Mouse FAS/Cd95 ELISA Kit (orb50156)

Description	Mouse FAS/Cd95 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q22-R169
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	7116 MW
Uniprot ID	Q8C350
Sensitivity	3 pg/ml
Expiration Date	Please enquire.



Mouse FAS
PicoKine
ELISA Kit
standard
CU...