


Product Datasheet

Human IL-5/Interleukin-5 ELISA Kit (orb50157)

Description	Human IL-5/Interleukin-5 ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: I20-S134
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	15238 MW
Uniprot ID	P05113
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 12.5 pg/ml.
Expiration Date	Please enquire.



Human IL-5
PicoKine
ELISA Kit
standard
C...