


Product Datasheet

Human KIM1 / TIM-1 ELISA Kit (orb50163)

Description	Human KIM1 / TIM-1 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S21-T288
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum, plasma (heparin, EDTA) and urine.
MW	38720 MW
Uniprot ID	Q96D42
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 31.2 pg/ml.
Expiration Date	Please enquire.



Human KIM1
PicoKine
ELISA Kit
standard c...