

Product Datasheet

phospho-Androgen Receptor Ser94 antibody (Biotin) (orb501639)

Description

phospho-Androgen Receptor Ser94

Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Androgen Receptor around the phosphorylation site of Ser94
Target	phospho-Androgen Receptor (Ser94)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	b3
Uniprot ID	P10275
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.