


Product Datasheet

Mouse MMP-3/Stromelysin-1 ELISA Kit (orb50168)

Description	Mouse MMP-3/Stromelysin-1 ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Y18-C477
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	53845 MW
Uniprot ID	P28862
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:40, concentration in serum and plasma is around 60 ng/ml.
Expiration Date	Please enquire.



Mouse
MMP-3
PicoKine
ELISA Kit
standard ...