


Product Datasheet

Mouse KIM1 / TIM1 ELISA Kit (orb50171)

Description	Mouse KIM1 / TIM1 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Y22-T212
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	30968 MW
Uniprot ID	Q3V033
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 0-208 pg/ml.
Expiration Date	Please enquire.



Mouse KIM1
PicoKine
ELISA Kit
standard c...