


Product Datasheet

Human BMP-5 ELISA Kit (orb50174)

Description	Human BMP-5 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A317-H454
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and bone tissue.
MW	51737 MW
Uniprot ID	P22003
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 156 pg/ml.
Expiration Date	Please enquire.



Human BMP-5 PicoKine ELISA Kit standard ...