

Product Datasheet

Insulin Receptor (phospho-Tyr999) antibody (Biotin) (orb501743)

Description

Insulin Receptor (phospho-Tyr999)

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Biotin
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Insulin Receptor around the phosphorylation site of Tyr999
Target	Phospho-Insulin Receptor (Tyr999)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	3
Uniprot ID	P06213
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500
Expiration Date	12 months from date of receipt.