


Product Datasheet

Human IL-1RA/IL1RN ELISA Kit (orb50177)

Description	Human IL-1RA/IL1RN ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: R26-E177
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	20055 MW
Uniprot ID	P18510
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 800 pg/ml.
Expiration Date	Please enquire.



Human IL-1RA
PicoKine
ELISA Kit
standard...