


Product Datasheet

Mouse MCP-1 ELISA Kit (orb50178)

Description	Mouse MCP-1 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q24-R96
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	16326 MW
Uniprot ID	P10148
Sensitivity	6 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum is 51-259 pg/ml.
Expiration Date	Please enquire.



Mouse MCP-1 PicoKine ELISA Kit standard ...