


Product Datasheet

Mouse BAFF ELISA Kit (orb50179)

Description	Mouse BAFF ELISA Kit
Reactivity	Mouse
Range	46.9 pg/ml - 3,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A127-L309
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	34192 MW
Uniprot ID	Q3KP92
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:50, concentration in serum and plasma is 20-30 ng/ml.
Expiration Date	Please enquire.



Mouse BAFF
PicoKine
ELISA Kit
standard c...