


Product Datasheet

Human Angiogenin/ANG ELISA Kit (orb50181)

Description	Human Angiogenin/ANG ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q25-P147
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	16550 MW
Uniprot ID	P03950
Sensitivity	12 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is 156-524 ng/ml.
Expiration Date	Please enquire.



Human ANG
PicoKine
ELISA Kit
standard
CU...