


Product Datasheet

Human EG-VEGF/PK1/PROK1 ELISA Kit (orb50185)

Description	Human EG-VEGF/PK1/PROK1 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A20-F105
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	11715 MW
Uniprot ID	P58294
Sensitivity	7 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 100 pg/ml.
Expiration Date	Please enquire.



Human EG-VEGF
PicoKine
ELISA Kit
standar...