


Product Datasheet

Mouse Resistin ELISA Kit (orb50186)

Description	Mouse Resistin ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S21-S114
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	12492 MW
Uniprot ID	Q5BMX4
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:100, concentration in serum and plasma is 14-100 ng/ml.
Expiration Date	Please enquire.



Mouse
Resistin
PicoKine
ELISA Kit
standa...