


Product Datasheet

Human CD40/TNFRSF5 ELISA Kit (orb50191)

Description	Human CD40/TNFRSF5 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E21-R193
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, urine and human milk.
MW	30619 MW
Uniprot ID	P25942
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 90 pg/ml.
Expiration Date	Please enquire.



Human
CD40/TNFRSF5
PicoKine ELISA
Kit st...