


Product Datasheet

Mouse CD40/TNFRSF5 ELISA Kit (orb50192)

Description	Mouse CD40/TNFRSF5 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V24-R193
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	32093 MW
Uniprot ID	P27512
Sensitivity	1 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 39-183g/ml.
Expiration Date	Please enquire.



Mouse
CD40/TNFRSF5
PicoKine ELISA
Kit st...