


Product Datasheet

Human EMMPRIN/CD147/BSG ELISA Kit (orb50194)

Description	Human EMMPRIN/CD147/BSG ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T25-H205
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, urine and human milk.
MW	42200 MW
Uniprot ID	P35613
Sensitivity	2 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 40 pg/ml.
Expiration Date	Please enquire.



Human
Emmprin
PicoKine
ELISA Kit
standar...